



U.S. Micro-Solutions, Inc. * 302 Unity Plaza * Latrobe, PA 15650
 Phone: (724) 853-4047 Fax: (724) 853-4049 AIHA-LAP, LLC EMLAP # 103009
www.usmslab.com

Customer Name: U.S. Micro-Solutions, Inc. Sample Date: September 10, 2020
 Customer Address: 302 Unity Plaza Date Received: September 11, 2020
 Latrobe, PA 15650 Date of Report: September 14, 2020
 Customer Phone: (724) 853-4047 Fax: (724) 853-4049
 PO Number: Attention:
 Project Name/Number: Swab Sample Report

Customer sample numbers below are uniquely identified by prefixing Laboratory # 12345-20

24 Hour Culturable Swab Sample(s) (Bacteria) - Analytical Method MIC 03 MIC 04 MIC 18

Sample Number	Media	Sample Description	Results of Microbial Analysis
1	TSA MAC	A	Total Bacterial Count < 100 CFU/swab

AS (TSA) = 100 CFU/swab No growth
 AS (MAC) = 2 CFU/swab

Total Raw Count: <1
 Analytical Sensitivity: 100 CFU/swab

2	TSA MAC	B	Total Bacterial Count 1,000 CFU/swab
---	------------	---	--------------------------------------

AS (TSA) = 100 CFU/swab Gram-negative bacteria isolated
 AS (MAC) = 2 CFU/swab

Sample Report

Total Raw Count: 10
 Analytical Sensitivity: 100 CFU/swab

3	TSA MAC	C	Total Bacterial Count 100 CFU/swab
---	------------	---	------------------------------------

AS (TSA) = 100 CFU/swab No growth of gram-negative bacteria
 AS (MAC) = 2 CFU/swab

Total Raw Count: 1
 Analytical Sensitivity: 100 CFU/swab

Samples are in good condition unless otherwise noted. Results relate only to the samples tested as received. Results are reported as calculated. For biological data, the first and/or second digit should be considered significant.

When providing duplicates of this report, the document should be provided in total and not in section in accordance with AIHA-LAP, LLC. Any unauthorized or improper disclosure, copying, distribution, use, or falsification of these results is prohibited. USMS shall have no liability to the Customer or the Customer's customer for opinions stated, recommendations made, actions taken, or conduct implemented based on the test results reported.

Technical Manager:

Deanna L. Kiska

Deanna L. Kiska, Ph.D.