

Enverify™ Technician Sampling Performance Report

ABC Pharmaceutical Company
 2107 Mulberry Street
 Yankton, SD 57078
 Phone #: 605-555-0178

Sampling Date: 11/28/2022
 Test Date: 11/30/2022
 Report Date: 12/2/2022

Technician: Liam P.

Sampling Device	Manufacturer	Cat#	Lot#	Expiration Date
Contact plate	Sampler Co.	KLZ-387	KP94346	6/3/2024

Enverify™ Test Surface	Lot#	Average CFU/Live Sample
	E25922C-220907	60

Results:

Table 1. Live Sample Results.

Sample Description	CFU/Sample	% Recovery
Live Sample 1	34	56.7
Live Sample 2	28	46.7
Live Sample 3	28	46.7
Mean	30	50.0
Standard Deviation	3.5	5.8

Table 2. Blank Sample Results.

Sample Description	CFU/Sample
Blank Sample 1	0
Blank Sample 2	0

Table 3. Summary of Results.

Test	Result	Specification	PASS/FAIL
Live Sample % Recovery (Mean)	50.0%	≥ 20.0%	PASS
Blank (Mean)	0 CFU	0 CFU	PASS