



## Enverify™ Result Notification – Pass

Please find the attached performance report for your Enverify™ Sampling Competency project. Your digital certificate and badge will be sent separately from [certificates@stratixlabs.com](mailto:certificates@stratixlabs.com)

## Enverify™ Technician Sampling Performance Report

ABC Pharmacy  
 2107 Mulberry Street  
 Yankton, SD 57078  
 Phone #: 605-555-0178

Sampling Date: 10/17/2024  
 Test Date: 10/21/2024  
 Report Date: 10/21/2024  
 Report Reference: 10-03790-24

**Participant: Liam P.**

Sampling Device	Manufacturer	Cat#	Lot#	Expiration Date
Contact plate	bioMérieux	N/A	1010722200	12/22/2024

Enverify™ Test Surface	Lot#	Average CFU/Live Sample
	E25922C-20240111	71.1

**Results:**

**Table 1. Live Sample Results.**

Sample Description	CFU/Sample	% Recovery
Live Sample 1	40	56.3
Live Sample 2	37	52.1
Live Sample 3	32	45
Mean	36.3	51.1
Standard Deviation	4.0	N/A

**Table 2. Blank Sample Results.**

Sample Description	CFU/Sample
Blank Sample 1	0
Blank Sample 2	0

**Table 3. Summary of Results.**

Test	Result	Specification	PASS/FAIL
Live Sample % Recovery (Mean)	51.1%	≥ 20.0%	PASS
Blank (Mean)	0 CFU	0 CFU	PASS